









Instructions for Use

FORMFIT[®] POST-OP KNEE

		3
EN	Instructions for Use	4
DE	Gebrauchsanweisung	7
FR	Notice d'utilisation	10
ES	Instrucciones para el uso	13
IT	Istruzioni per l'uso	16
NO	Bruksanvisning	19
DA	Brugsanvisning	22
SV	Bruksanvisning	25
EL	Οδηγίες Χρήσης	28
FI	Käyttöohjeet	31
NL	Gebruiksaanwijzing	34
PT	Instruções de Utilização	37
PL	Instrukcja użytkownika	40
CS	Návod k použití	43
TR	Kullanım Talimatları	46
RU	Инструкция по использованию	49
JA	取扱説明書	52
ZH	中文说明书	55
KO	사용 설명서	59
SK	Návod na použitie	62
ET	Kasutusjuhend	65
RO	Indicații pentru Utilizare	68
AR	تعليمات الاستخدام	71

ENGLISH

PRODUCT DESCRIPTION

The Formfit® Post-Op Knee 16" brace controls the range of motion (ROM) or immobilizes the knee after surgery or injury (**Figure 1**). Two sizes are available, regular and XL. Replacement softgoods and accessories can be ordered separately.

Sizes

Type	Thigh Circumference	Order Number
Formfit® Post-Op Knee 16" Reg	Up to 27" / 69 cm	FPOK403300
Formfit® Post-Op Knee 16" XL	27" – 35" / 69–89 cm	FPOK403301

Replacement Softgoods

Type	Order Number
Regular	403275
XL	403277

Accessories

Type	Order Number
ROM stop locks (set of 20 lock clips)	B-8380005

INDICATIONS FOR USE

Post-operative/post-injury ROM control or immobilization.

Ligament and fracture-related surgeries and when controlled range of motion is necessary.

CONTRAINDICATIONS

None known.

INTENDED USE

The Formfit® Post-Op Knee is intended to be used for knee immobilization and range of motion restriction.

GENERAL SAFETY INSTRUCTIONS

WARNING: Do not use for multiple patients. Only use the device for a single patient.

CAUTION: If you experience increased pain, swelling, sensation changes or any adverse reactions while using this product, immediately consult your medical professional.

USA Only: Federal law restricts this device to sale by or on the order of a licensed health care practitioner.

USING THE BRACE

Size Selection

Measure the thigh circumference 6" / 15 cm above the mid patella. Select a size that agrees with the measurement. Refer to the "Sizes" chart in the section "Product Description" to find the correct size.

Donning the Brace

WARNING: Make sure to keep the leg straight when you slide the foam pads and straps under it.

1. Unfasten the leg straps.
2. Fully pull the open ends of the straps through the double D rings to remove them from the left strut.
3. Remove the foam pads from the two struts.
4. Fully open the foam pads.
5. Carefully slide the foam pads under the thigh and calf (**Figure 2**).
6. Fasten the foam pads around the leg (**Figure 3**).
7. Put the struts on each side of the leg.
8. Align the top of the hinge pivot with the top of the patella (**Figure 4**).
9. Align the foam pads with the hook and loop fasteners on the inner surface of the struts.
10. Attach the struts to the foam pads on each side of the leg.
11. Carefully slide the straps under the foam wraps.
12. Put the straps back through the double D rings.
13. Pull and fasten the straps around the foam pads in the correct sequence (**Figure 5**).

CAUTION: Do not tighten the straps too much. If the straps are too tight, the blood circulation can be cut off.

14. If necessary, use bending tools to adapt the brace to the patient's leg shape (**Figure 6**).

Setting Range Of Motion

Extension can be limited at: -10°, 0°, 20°, 30°, 40°, 50°, 60°, 70°, 80°, 90°.

Flexion can be limited at: 10°, 20°, 30°, 40°, 50°, 60°, 70°, 80°, 90°, 100°, 110°, 120°.

Adjust ROM as follows (**Figure 7**):

Extension: Pull the black extension button up and set the angle.

Flexion: Pull the white flexion button up and set the angle.

Release a button to let it click into place. Put the stop locks onto the buttons to secure them.

Immobilizing the Knee

Use the drop lock to fully immobilize the knee. It can be set at -10°, 0°, 10°, 20°, 30°. An arrow shows the angle. The angle is shown by the scale on the strut, below the hinge.

Adjust immobilization as follows (**Figure 8**):

1. Push the blue drop-lock button up to unlock the hinge.
2. Move the brace until the arrow and the correct angle agree.
3. Push the blue drop-lock button down to lock the hinge angle.

Doffing the Brace

1. Unfasten the straps and the foam pads.
2. Remove the brace (**Figure 9**). You can put it back on again fully assembled.

PRODUCT CARE

The foam pads can be hand washed with mild soap. Leave to dry. Do not tumble dry.

PRODUCT MATERIAL LIST

Polyamide, polyamide with glass fiber, polyurethane, steel, stainless steel, aluminum, zinc.

DISPOSAL

All components of the product and packaging should be disposed of in accordance with respective national environmental regulations. Users should contact their local governmental office for information on how these items can be disposed of in an environmentally sound manner.

WARRANTY

Plastic parts and hinges 3 months warranty. Foam pads and leg straps 30 days. 6 months free from defects.

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