

REBOUND® AIR WALKER

Rebound Air Walker is a circumferential walker with fully integrated air pump and release valve that provides easy, tool-free inflation and deflation for fine-tuned, consistent compression. It features a host of details designed to provide added protection, simplified application and usage.



INDICATIONS FOR USE

Foot and ankle conditions that may benefit from motion restriction, such as:

- Grade 2 and 3 ankle sprains
- Stable fractures
- Tendon or ligamentous sprain
- During postoperative or rehab phase. If used with the Achilles Wedge Kit, Achilles ruptures can be treated as well.
- Rebound Air Walker Low Top
- Foot and ankle conditions that may benefit from motion restriction, such as:
- Stable foot fractures
- Tendon or ligamentous sprain
- During postoperative or rehab phase.

FEATURES AND BENEFITS

- Linked medial and lateral air bladders with a single inflation point provide balanced compression and are shaped to cradle the Achilles, holding the heel in place
- Rocker bottom promotes a natural, stable gait, while a trimmable, semi-soft toe bumper provides added protection
- Moisture-wicking softgood helps keep skin dry and comfortable
- Straps fasten in front for easy application, adjustment and removal
- EVA material provides greater shock absorption in the heel
- Fully-integrated air pump and release valve provides easy, tool-free inflation and deflation for fine-tuned, consistent compression
- Flex Edge™ overmolding conforms to each patient's unique anatomy to help limit the formation of pressure points
- Ventilated panels wick moisture to prevent microbial buildup and improve patient comfort

REQUIRED MEASUREMENTS

- Shoe size



REBOUND® AIR WALKER

ORDERING INFORMATION

REBOUND AIR WALKER

Size	Shoe Size UK	Shoe Size EU	Standard Liner / Standard	Standard Liner / Low Top	Slim Liner / Standard	Slim Liner / Low Top
X-Small	1.5 - 3.5	33 - 35.5	B-242900001E	B-242900061E	B-242900001SE	B-242900061SE
Small	4 - 6.5	36 - 39	B-242900002E	B-242900062E	B-242900002SE	B-242900062SE
Medium	7 - 10	39.5 - 44	B-242900003E	B-242900063E	B-242900003SE	B-242900063SE
Large	10.5 - 12	44.5 - 46.5	B-242900004E	B-242900064E	B-242900004SE	B-242900064SE
X-Large	12.5 - 14.5	47 - 50	B-242900005E	B-242900065E	B-242900005SE	B-242900065SE

ACCESSORIES

REPLACEMENT SOFTGOODS

Size	Standard Liner / Standard	Standard Liner / Low Top	Slim Liner / Standard	Slim Liner / Low Top
X-Small	B-418000101	B-418000161	B-418030101	B-418030161
Small	B-418000102	B-418000162	B-418030102	B-418030162
Medium	B-418000103	B-418000163	B-418030103	B-418030163
Large	B-418000104	B-418000164	B-418030104	B-418030164
X-Large	B-418000105	B-418000165	B-418030105	B-418030165

ACHILLES HEEL WEDGE KIT

Size	Description	Part#
Small	Wedge kit for walkers XS-S	AWK0001
Medium	Wedge Kit for walkers M	AWK0005
Large	Wedge kit for walkers L-XL	AWK0009
N/A	5 layer Wedge kit	AWK00015

EVENUP

Size	Part#
Small	10350
Medium	10351
Large	10351

WALKER COVERS

Size	Hygiene covers (2 pcs)	Durable weather cover	Disposable weather cover (20 pcs)
Small	B-418500602	B-418500502	B-418500512
Medium	B-418500603	B-418500503	
Large	B-418500604	B-418500504	B-418500514
X-Large	B-418500605	B-418500505	

OTHER ACCESSORIES

Description	Part#
Heel Wedge	B-242902511
Adhesive Insole for Wedge	WHW-W0107
Strap Extension	B-502100101
Rebound Hightop Sock SM/MD SM/MD	B-418500703
Rebound Hightop Sock LG/XL LG/XL	B-418500704
Rebound Lowtop Sock SM/MD SM/MD	B-418500903
Rebound Lowtop Sock LG/XL LG/XL	B-418500904