

# FORMFIT POST-OP KNEE

## INDICATIONS FOR USE

- Post-op/post-injury ROM control or immobilisation
- Indicated for ligament and fracture-related surgeries and when controlled range of motion is necessary

## FEATURES AND BENEFITS

- Breakaway struts for easy, tool free length adjustment
- OneTouch™ easy to adjust ROM hinge with lock-out capability from -10° to 30° (10° incremental)
- Extension control from -10° to 90° (10° incremental)
- Flexion control from 0° to 120° (10° incremental)
- Available in cool or foam models, 3 lengths, and two sizes (Regular and XL)



## REQUIRED MEASUREMENTS



Measure thigh circumference 15 cm (6") above mid-patella. Leg in extension.

Size	Thigh circumference 15 cm (6") above mid patella	
	cm	inch
XL	69 - 89	27 1/8 - 35
Regular	0 - 69	0 - 27 1/8

## ORDERING INFORMATION

### FULL FOAM, 61 CM (24")

Size	Part#
XL	FPOK403202
Regular	FPOK403200

### FULL FOAM, 61 CM (24"), BAG

Size	Part#
XL	FPOK403202I
Regular	FPOK403200I

### COOL FOAM, 61 CM (24")

Size	Part#
XL	FPOK403102
Regular	FPOK403100

### COOL FOAM, 61 CM (24"), BAG

Size	Part#
XL	FPOK403102I
Regular	FPOK403100I